BODY

COMFORT

SINCE 1992

CZ

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SK

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www.bodycomfort.biz

tel. č.: 0903 630 681

GB Instruction Manual + Guarantee sheet

Congratulation for choosing this paraffin wax spa BODY COMFORT BC 01, which has a high heat capacity to absorb and retain a great amount of heat. Paraffin wax spa takes advantage of paraffin heat retaining capacity in a simple process which is called heat transfer. Heat is transferred to the core of the affected area as paraffin area goes through a phase change. At the same time, the heat from paraffin opens pores and increases circulation in the skin, leaving skin softened and more radiant.

READ ALL INSTRUCTIONS in this manual before using and follow all warning instructions!

Functions + features of paraffin wax spa BODY COMFORT BC 01:

- soften and moisturize dry hands, elbows, feet and any other part of the body while putting a wax
- can ease pain in sore joints
- relief from tired muscles
- good by rheuma diseases

This handy paraffin wax spa turns your home into a „beauty salon+health spa“ .

Produkt description:

1. transparent plastic cover
2. inside of the wax tub
3. protective plastic grid
4. control panel
5. lower part of the wax heater
6. socket for electrical wire in the wax heater to plug in
7. power switch on the wax heater
8. button „II“ for disconnecting heating function and holding reached temperature only
9. start button „I“ for starting of heating function
10. left side power indicator light
11. right side heating indicator light
12. button selector for temperature increasing
13. button selector for temperature decreasing
14. display indicating temperature of the heating element in degree Celsius
15. plugged in electrical cable in the wax heater

This electrical appliance can be used only for the purpose described in this instruction manual. No other use is allowed !

When using this electrical appliance, especially when children or disable persons are present, all safety precautions must be always followed.

How to set up the paraffin wax spa:

- remove all the contents from the packaging. Open the bags of the paraffin wax spa and place the wax inside the tub.

Note: Make sure that the tub is clean and empty before putting the paraffin wax into the tub. This model of paraffin wax spa is designed for up to max. 5 kilograms of paraffin wax while you have to put inside minimum of 1 kilogram of wax to cover the bottom of the tub.

- This paraffin wax heater can be used/heated only with paraffin wax. No use of water or any liquid is allowed !

- plug in one side of the power cord into the socket in the wax heater /15/ and after that another side of the power cord to power socket on the wall, but be careful while this can be plugged only into the 230 Volt/ 50 Hz electrical net. Then turn the power switch on to the position „I“
- paraffin wax heater must be placed on the horizontal flat, fixed, heat-resistant surface
- for safety purpose you can use the protective plastic grid, but after the use this must be pulled out, while making the melting time much more longer
- for making on the power, press the power switch to the position „I“ and the right indicating power light will appear and on the display will appear the last adjusted temperature.
- for adjusting/changing the temperature of heating first press the start button „I“ on the controller and press up or down for selection of the temperature. Maximum temperature level in this wax heater is 85 degree Celsius. During increasing the temperature the heating indicating light is flashing. After reaching the selected temperature the light will stop flashing and will appear.
- for changing the selected temperature is always necessary to press the left up button „I“ followed by selecting the temperature and after selection of temperature once more pressing button „I“ on the controller.

- !!! Important notice !: Wax heater can be heated only when there is any paraffin wax inside ! Never heat up the wax heater without wax ! It is not allowed to put inside water or any liquid ! Do not heat up empty wax heater !
 - after reaching the adjusted temperature of the heating element/paraffin wax, the red indicating heating light will stop flashing, and will light only
 - paraffin wax can be used only melted, but not immediately after melting ! At that time it can be too hot, what can be dangerous while burning the skin ! Maximum temperature level of wax can be up to 60 degree Celsius and less ! Otherwise the skin can be burned !
 - fastest melting time during the first melting of 1 kilogram low melting point wax „46“ inside the tub is by adjusting this wax heater on the maximum of 85 degree Celsius heating, around 50 minutes. By following next melting, the time is up to 55 minutes.
 - melting time on average depends on many different factors, beside others from amount of the wax, melting point of the wax and maximum heating temperature. For example 2 kilograms of low melting point „46“ paraffin wax BODY COMFORT melting time is only up to-around 60 minutes !
 - !!! Important notice !: do not immerse hands, feet or any part of the body or skin inside the wax heater immediately after melting the wax ! At that time the wax temperature will be too high and can burn the skin ! The wax temperature at that time could have also around 85 degree Celsius what's very dangerous for the skin !
 - Paraffin wax can be used on the skin only after lowering the temperature of the wax on the accepted temperature of less than 60 degree Celsius, optimal would be around 50-55 degree Celsius. This can be reached by setting-adjusting the temperature on the controller to the level of around 55 degree Celsius and waiting some few time – usually around 25-35 minutes, until the wax temperature drops down. Before using is suggested to prove the wax temperature by checking it with the finger
 - purpose of the protective plastic grid is burning protection while putting the hands, elbows, arms inside of the tub
 - protective plastic tube must be pulled out from the tub after using-therapy
 - temperature of the paraffin wax by this wax heater can be adjusted just after the wish of the user up to 85 degree Celsius, only is needed to be considered, that the temperature over 60 degree Celsius is dangerous for human skin and can burn the skin !
 - paraffin wax included in this wax heater is a specially approved cosmetic+medical wax BODY COMFORT with low melting point of around 46 degree Celsius, means to use is melted the temperature controller needs to be adjusted to around 50-55 degree Celsius
 - function of „II“ button on the controller is immediate stop of heating and holding the warmth of the wax at the reached level
 - ! Important notice: when heating will reach the adjusted level, heating will automatically stop and will be holding the range of the adjusted preselected temperature. It means, it will not stay exactly on the settled level, but will go some few degree up for warming and after that will fall down few degree also, after which will be automatically heated repeatedly to, that it will move in the temperature level around +/- max. 5 degree over/under the adjusted level, what is usual
 - for the safe use the hands, elbows, arms + others can be putted inside of the wax tube only dry, clean, without any metals + others like rings etc. and only after some 3-5 minutes after displayed maximum 60 degree Celsius on the display, or less
 - for your safety please try your best to avoid to touch the aluminium surface of the tub, whether on the bottom or on the side walls, to prevent eventually burning of the skin
 - this paraffin wax heater BODY COMFORT BC 01 is good only for human use, individual or also professional use !
- Important notice:

It's strongly recommended the testing of the paraffin wax temperature on the inside of the wrist area before dipping on your hand, elbow or foot entirely !

Place the grate into the unit only when the paraffin wax is melted.

How to use the paraffin wax treatment:

First you have to pull out all the jewellery + bracelets + rings + diamonds + watches + others from the area which will be treated. Clean the area of the body part and skin to be treated with warm water and soap, and then apply a light coat of lotion on the skin before dipping into the wax. It will allow the wax to be removed easier and help moisturize the skin.

Important notice: you have to test the paraffin wax temperature on the inside of your wrist to make sure the temperature of the wax is comfortable before immersing your hand, elbow, foot or any other part of your body, into the paraffin wax.

Immerse the hand or elbow into the paraffin wax, hold it there for some few seconds only and then remove it out. Before redipping you have to wait some 5 seconds /until the paraffin wax no longer appears shiny/. The point is to dip and pull out of the wax the body part to build up layers of paraffin wax.

Important notice: Do not leave any part of your body or skin in the tub for long time, do not soak, dip quickly in and out !

Repeat this process 3 to 5 times ! For relief of joint and stiffness dip 10 to 15 times to form a very thick glove. For very dry cracked skin apply at least 5 layers of paraffin.

Cover the treated area with a protective plastic liner and then slip on a mitt/boot. This retains heat for a longer time treatment. Allow the paraffin wax to remain in the place until the heat goes away = about 15 to 20 minutes.

When ready to remove the paraffin wax, first remove the mitts and plastic liners, then slide off the paraffin wax. Discard used plastic liners.

Normally discard used wax. But when you are the only person using this paraffin wax spa and the wax, then the wax could be reused, just put it back into the tub to remelt.

Important notice: you have to follow all the hygienical standards and regulation valid in your country to protect yourself and to avoid any health or safety problems!

Maintenance:

Unplug the electrical appliance and remove grate. Allow paraffin wax to solidify.

Push on one end of the paraffin wax spa block and it will tip up. Lift the block out of the machine or dispose of it. ! Do not use sharp subjects or objects to remove the paraffin wax block. This may scratch or damage the tub. If the paraffin wax does not remove easily, plug the machine in and allow to heat for about 15 minutes, then remove the paraffin wax block and unplug the appliance.

Wipe inside of the tub with the disposable cloth to remove residual dirt and rest of paraffin wax.

!!! IMPORTANT NOTICE: do not pour water into the paraffin wax heater or do NOT immerse this appliance into the water !!! You can use a damp soapy cloth to wipe clean this appliance.

This paraffin wax spa must be stored in a cool and dry place.

Technical data:

- power source electrical net 230 Volt/50 Hz
- Power consumption: max. : 300 Watt only
- CE and ROHs tested and certified

IMPORTANT SAFETY INSTRUCTIONS and PRECAUTIONS:

Please read all these instructions and follow them fully:

- To reduce the risk of electric shock:
- do not leave this appliance unattended while switched on
- do not place this appliance from where it could fall or to be pulled into the tub or sink. Do not place in or drop into water or any other liquid
- always before cleaning unplug this appliance from the electrical net-outlet
- do not plug or unplug this appliance when hands are wet.

WARNINGS:

To reduce the risks of burns, fire, electric shock or injury to persons you have to follow all these instructions:

- this appliance is no toy, can not be used by children, close supervision is necessary when this appliance is used by, on or near children, invalids or disabled persons
- this appliance can NOT be used by any person suffering from any physical ailment that would limit the users capacity to operate the controls
- use only attachments and wax provided with the unit or delivered by manufacturer of this paraffin wax spa
- do not use attachments NOT provided with this unit

- do NOT operate this appliance if it has a damaged cord or plug, or if it is not working properly, if it has been dropped into the water or any liquid
- keep electrical power cord away from heated or hot surfaces
- do not carry this appliance by power supply cord or do NOT use the power supply cord as a handle !
- do not use this product outdoors !This product is designed and developed for indoor use only
- do not let the power cord hang over the edge of the table or any counter
- do not place this product on or near any hot place, hot gas, flames, hot electric element, on, inside or near a hot or heated oven.Do NOT place this produkt on the top of any other appliance
- do not operate where aerosol spray products are being used or where oxygen is being administered
- this appliance is not intended for childrens use !
- never use this product while sleeping or dropsy
- never use this product while in bed
- do not apply the hot wax near eyes or other sensitive areas of the body
- do NOT touch the hot surfaces
- NEVER use any sharp instruments to remove solidified paraffin wax from the unit
- use appliance on a sturdy level surface only
- always store this appliance in a cool and dry place

SPECIAL WARNINGS:

- in case you have, had, or have had ANY skin or other health problems, diseases and or similar health problems you MUST consult your doctor-physician BEFORE first use of this appliance!
- use this appliance only in sitting position – while seated.Do NOT stand on or in the tub !
- do not apply paraffin wax on any open cuts or wounds, warts, inflamed skin, or when the sensation of any extremities is decreased
- you have to consult your doctor-physician, if you have any questions or concerns about acute inflammation and peripheral vascular disease or decreased sensation of any extremity before using
- discontinue using and consult your physician if the paraffin wax feels to be too hot or if it does not feel warm enough.This may indicate a possible physical problem
- immediately discontinue use of this appliance and consult your doctor-physician in case of dermatitis,discomfort, irritation, or pain occurs from continued paraffin treatments
- do NOT overfill the appliance with paraffin wax and always keep the paraffin wax level about less than the top allowed and marked inside of the tub
- NEVER use this appliance WITHOUT paraffin wax !!!
- do NOT heat paraffin wax in any of its packaging !! Paraffin wax can be heated ONLY without any packaging!
- paraffin wax can not be heated over an open flame, on a range,stove or in a microwave !
- do not pour the melted paraffin wax from the tub !
- do NOT use any accessories, wax or others, different from delivered by this manufacturer.Break of this rule can cause breaking of the guarantee while manufacturer of this produkt can NOT be responsible for destroying of this product while using unknown attachments, spare parts, different wax etc.It is **STRONGLY RECOMMENDED** to use all the manufacturers parts, attachments like wax, gloves etc.
- in ANY case of any problems, please contact your next local authorized service to solve your problem as soon as possible and in competent and qualified way.Any use of any uncompetent and unauthorized steps into this product through unauthorized staff will cancel any guarantee rights,
- the user of this appliance have to follow all the safety regulations valid in the country of use and is fully responsible for this

Guarantee sheet

Manufacturer of paraffin wax spa BODY COMFORT model BC 01 is fully responsible for all the production damages within 2 years from the date of the purchase by individual use and 6 months after the date of purchase in case of any professional use. During this date period manufacturer will solve all the problems of the product caused in production. This guarantee does NOT include any damages caused by user, unauthorized staff or similar, shipping cost to the service. Manufacturer is also not responsible for damages caused through use of this product for different purpose as described in this instruction manual. In any case of any NON authorized hit-intervention into this product, fall out any right or requirement for any guarantee. To apply sucesfully for a guarantee service you need to present the original bill of the authorized dealer/seller.

Any reclaim rights or requirements kindly please contact following authorized service: our authorized service staff is in Czech republic + Slovakia, adresses are in disposal on the previous pages.

Also in any case of any unclearancies concerning using this wax heater please do not hesitate to send us email to: info@bodycomfort.biz or eventually to call to the mentioned phone numbers.

Date of the purchase:

Place, stamp and signature of the authorized seller:

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Copyright © BODY COMFORT 08/ 2014 – 2024