

Cryospray SMT CS 1-SM



Type CS 1-SM

CE 1014

Czech Ministry of Health
Free Sale Certificate (FSC)
number MZDRX00GZ54Y

Cryosprays SMT CS 1-SM are designed for effective treatment of a wide range of skin lesions: Common verruca, seborrhic veruca, actinic keratosis, papilloma, lentiga and other common lesions, inclusive haemorrhoids.

This is the coldest (-196°C/-320F), the most effective, the most universal and accessible cryosurgical method. Easy and rapid. 3 to 30 seconds is sufficient for majority of operations. It is more effective and much more safe for patients than liquid nitrogen applied with swab. In contrast to scalpel surgery or electrosurgery: No or minimal pain, no blood, no offensive odour, usually no anesthesia, no bandages, no assistant needed, no medical waste. Excellent cosmetic effect. None or minimum scarring.

Cryospray SMT CS 1-SM je určen na efektivní léčení široké škály kožních lézí: Běžné veruky, seborhoické veruky, aktinické keratomy, papilomy, lentiga, basaliomy a další běžné kožní léze. Kromě toho lze s výhodou použít i na chirurgické zákroky, například na hemeroidy. Jedná se o nejchladnější -196°C, a tedy o nejúčinnější a nejuniverzálnější dostupnou kryochirurgickou metodu. Aplikace je jednoduchá a rychlá - pro většinu zákroků postačí 3 - 30 sekund. Mnohem účinnější než dusík aplikovaný „štetěčkami“. Na rozdíl od chirurgie skalpelem nebo elektrochirurgie: Žádná nebo minimální bolest, žádný zápach, nekrvavý zákrok, není třeba obvaz nebo bandáž, většinou není třeba anestezie, není třeba sestry pro asistenci, žádný medicínský odpad, výborný kosmetický efekt, minimální nebo žádná jizva.

TECHNICAL PARAMETERS

Coolant	liquid nitrogen (- 196°C)
Capacity	0,2l
Weight without LN2	325 g
Weight with LN2	460 g
Height	175 mm
Diameter	58 mm
Pressure	75 kPa
Insulation	vacuum

TECHNICKÉ PARAMETRY

Chladičí médium	kapalný dusík (- 196°C)
Objem nádoby	0,2 l
Hmotnost bez náplně	325 g
Hmotnost s náplní	460 g
Výška	175 mm
Průměr	cca 58 mm
Pracovní tlak	cca 75 kPa
Izolace	vakuová